



FELINE SERUM ALBUMIN
STANDARD GRADE POWDER

TDS- FSA62

Effective Date: February 19,2019

Rev. 000

Page 1 of 1

METHOD OF PURIFICATION: Prepared from feline serum by heat-shock method.

SPECIFICATIONS:

Form:	Powder
Color:	Tan to amber
Source:	Feline serum
Moisture:	$\leq 5.0\%$
Product pH (4% solution in DI H₂O):	7.0 ± 0.2
Albumin - Purity: (Electrophoresis)	$\geq 96\%$
Protein level:	$\geq 3.6 \text{ g/dl } (\geq 96\%)$

CATALOG NUMBERS:

FSA62-0050

FSA62-0100

SIZES:

50 mg

100 mg

SHIPPING CONDITIONS: Ambient temperature

STABILITY AND STORAGE:

Store lyophilized powder refrigerated at $5^{\circ}\text{C} \pm 3^{\circ}\text{C}$ in sealed containers. Store reconstituted material frozen at $-20^{\circ}\text{C} \pm 5^{\circ}\text{C}$. Avoid repeated freezing and thawing.

The material for this product was processed in the USA from the blood of healthy animals that had not been inoculated with or exposed to any infectious agents. During processing this material was subjected to a pH of $\leq 5^{\circ}\text{C}$ and $\geq 65^{\circ}\text{C}$ for at least 3 hours. At no time during collection or processing was the material commingled with any other material of animal origin.

This product is not to be used for direct or indirect exposure or inoculation into laboratory and domestic animals.

This product is sold for use in IVD, research, further manufacturing, or technical purposes only, and is not intended for direct use in diagnostic or therapeutic applications.