



**BOVINE SERUM  
DELIPIDIZED  
STERILE FILTERED**

TDS- SBUDP30

Effective Date: March 13, 2025

Rev. 000

Page 1 of 1

**SPECIFICATIONS:**

Source:	Bovine Serum
Hemoglobin:	<40 mg/dl
Protein:	≥5.0 g/dl
pH:	6.5 – 8.5
Osmolality:	260 – 370 mOsm/Kg
Filtration:	0.2 micron (Absolute)
Sterility:	No Growth

**CATALOG NUMBER:**

SBUDP30-0100  
SBUDP32-0100  
SBUDP30-0500  
SBUDP32-0500

**SIZES:**

100ml  
100ml Heat Inactivated  
500ml  
500ml Heat Inactivated

**SHIPPING CONDITIONS:** <0°C

**STABILITY AND STORAGE:**

This serum is stable for ≥ 5 years when stored -20° C ± 5°. Avoid repeated freezing and thawing.

The raw materials for these products were collected from cattle 30 months and younger in the USA in abattoirs registered with the USDA. The animals received ante and post-mortem inspections under a veterinarian's supervision and were apparently free from infectious and contagious diseases. At no time during collection or processing was the material commingled with any other material of animal origin.

This product is sold for use in vitro diagnostic, research, further manufacturing, or technical purposes only, and is not intended for direct use in diagnostic or therapeutic applications.