



**CANINE SERUM
STERILE FILTERED**

TDS-SCA30

Effective Date: May 22, 2019

Rev. 000

Page 1 of 1

SPECIFICATIONS:

Source:	Dog serum, mixed sex, mixed breed
Hemoglobin:	<30.00 mg/dl
Protein:	>5.0 g/dl
pH:	6.8 - 8.2
Osmolality:	240 – 360 mOsm/Kg
Filtration:	0.2 micron (Absolute)
Sterility:	No Growth

CATALOG NUMBER:

SCA30-0050
SCA32-0050
SCA30-0100
SCA32-0100
SCA30-0500
SCA32-0500

SIZES:

50ml
50ml Heat Inactivated
100ml
100ml Heat Inactivated
500ml
500ml Heat Inactivated

SHIPPING CONDITIONS: <0°C

STABILITY AND STORAGE:

This serum is stable for ≥ 3 years when stored $-20^{\circ} \text{C} \pm 5^{\circ}$. Avoid repeated freezing and thawing.

The material for this product was processed in the USA from the blood of healthy dogs that had not been inoculated with or exposed to any infectious agents. At no time during collection or processing was the material commingled with any other material of animal origin.

All of the dogs are vaccinated against rabies, canine parvovirus, canine distemper virus, and canine viral hepatitis virus. This product is not to be used for direct or indirect exposure or inoculation into laboratory and domestic animals.

This product is sold for use in IVD, research, further manufacturing, or technical purposes only, and is not intended for direct use in diagnostic or therapeutic applications