



**DONOR EQUINE SERUM  
DELIPIDIZED  
STERILE FILTERED**

**TDS-SEDP30**

Effective Date: May 22, 2019

Rev. 000

Page 1 of 1

**SPECIFICATIONS:**

<b>Source:</b>	Donor Equine Serum
<b>Hemoglobin:</b>	<10 mg/dl
<b>Protein:</b>	≥4.0 g/dl
<b>pH:</b>	7.5 - 8.9
<b>Osmolality:</b>	240 – 360 mOsm/Kg
<b>Filtration:</b>	0.2 micron (Absolute)
<b>Sterility:</b>	No Growth

**CATALOG NUMBER**

<b>SEDP30-0100</b>	<b>100ml</b>
<b>SEDP32-0100</b>	<b>100ml Heat Inactivated</b>
<b>SEDP30-0500</b>	<b>500ml</b>
<b>SEDP32-0500</b>	<b>500ml Heat Inactivated</b>

**SHIPPING CONDITIONS:** <0°C

**STABILITY AND STORAGE:**

This serum is stable for ≥ 5 years when stored -20° C ± 5°. Avoid repeated freezing and thawing.

The serum was produced from animals in the USA. This serum was collected in USDA approved collection facilities. The animals received inspections under a veterinarian's supervision and were apparently free from infectious and contagious diseases. At no time during collection or processing was the material commingled with any other material of animal origin.

This product is sold for use in IVD, research, further manufacturing, or technical purposes only, and is not intended for direct use in diagnostic or therapeutic applications.