

	<b>FETAL BOVINE SERUM STERILE FILTERED- US ORIGIN</b>	<b>TDS-SFBU30</b>
<b>Effective Date:</b>	<b>June 20,2019</b>	<b>Rev. 000</b>
		<b>Page 1 of 1</b>

**SPECIFICATIONS:**

<b>Source:</b>	Fetal Bovine Serum
<b>Hemoglobin:</b>	≤ 20 mg/dl
<b>Protein:</b>	≥3 g/dl
<b>pH:</b>	7.0-8.0
<b>Osmolality:</b>	280-340 mOsm/kg
<b>Endotoxin level:</b>	≤ 10 EU/ml
<b>Mycoplasma:</b>	Negative
<b>Cytotoxicity test:</b>	Negative
<b>Sterility:</b>	No growth

**CATALOG NUMBER:**

SFBU30-0010x10  
SFBU32-0010x10  
SFBU30-0050  
SFBU32-0050  
SFBU30-0100  
SFBU32-0100  
SFBU30-0500  
SFBU32-0500

**SIZE:**

10x10ml  
10x10ml Heat Inactivated  
50ml  
50ml Heat Inactivated  
100ml  
100ml Heat Inactivated  
500ml  
500ml Heat Inactivated

**SHIPPING CONDITIONS:** <0°C

**STABILITY AND STORAGE:**

For long term storage, store frozen at -20° ± 5°C. Avoid repeated freezing and thawing.

The raw materials for these products were collected from animals in the USA in abattoirs registered with the USDA. The animals received ante and post-mortem inspections under a veterinarian's supervision and were apparently free from infectious and contagious diseases. At no time during collection or processing was the material commingled with any other material of animal origin.

This product is sold for use in vitro diagnostic, research, further manufacturing, or technical purposes only, and is not intended for direct use in diagnostic or therapeutic applications.