



**DONOR GOAT SERUM  
DELIPIDIZED  
STERILE FILTERED**

**TDS-SGDP30**

Effective Date: March 13, 2025

Rev. 000

Page 1 of 1

**SPECIFICATIONS:**

Source:	Donor Goat Serum
Hemoglobin:	<40 mg/dl
Protein:	≥5.0 g/dl
pH:	6.5 – 8.5
Osmolality:	260 – 370 mOsm/Kg
Filtration:	0.2 micron (Absolute)
Sterility:	No Growth

**CATALOG NUMBER:**

SGDP30-0100  
SGDP32-0100  
SGDP30-0500  
SGDP32-0500

**SIZES:**

100ml  
100ml Heat Inactivated  
500ml  
500ml Heat Inactivated

**SHIPPING CONDITIONS:** <0°C

**STABILITY AND STORAGE:**

This serum is stable for ≥ 3 years when stored -20° C ± 5°. Avoid repeated freezing and thawing.

The serum was produced from animals in the USA. This serum was collected in USDA approved collection facilities. The animals received inspections under a veterinarian's supervision and were apparently free from infectious and contagious diseases. At no time during collection or processing was the material commingled with any other material of animal origin.

This product is sold for use in vitro diagnostic, research, further manufacturing, or technical purposes only, and is not intended for direct use in diagnostic or therapeutic applications.