



**BOVINE NEWBORN CALF SERUM  
STERILE FILTERED, HEAT  
INACTIVATED**

**TDS-SNBU32**

Effective Date: 1/10/2025

Rev. 000

Page 1 of 1

**SPECIFICATIONS:**

<b>Source:</b>	Bovine Newborn Calf Serum
<b>Hemoglobin:</b>	<30 mg/dl
<b>Protein:</b>	>3.5 g/dl
<b>pH:</b>	6.8 - 8.2
<b>Osmolality:</b>	260 – 360 mOsm/Kg
<b>Filtration:</b>	0.2 micron (Absolute)
<b>Sterility:</b>	No Growth

**CATALOG NUMBER:**

**SNBU32-0100**

**SNBU32-0500**

**SIZES:**

**100ml**

**500ml**

**SHIPPING CONDITIONS:** <0°C

**STABILITY AND STORAGE:** This serum is stable for  $\geq 5$  years when stored  $-20^{\circ} \text{C} \pm 5^{\circ}$ . Avoid repeated freezing and thawing.

The raw materials for these products were collected from calves 10 days and younger in the USA in abattoirs registered with the USDA. The animals received ante and post-mortem inspections under a veterinarian's supervision and were apparently free from infectious and contagious diseases. At no time during collection or processing was the material commingled with any other material of animal origin.

This product is sold for use in vitro diagnostic, research, further manufacturing, or technical purposes only, and is not intended for direct use in diagnostic or therapeutic applications.